



Therapeutic Oligonucleotides

Agile Support fom Discovery to Preclinical Success



A wide range of synthetic ssRNAs are available modified and unmodified.



Specifications

- **Length:** From 5 to 200 bases

- **Synthesis Scale:** mg to gram quantities

- **Backbone:** RNA, PS, PMO, 2'O-Me RNA, 2'O-MOE RNA, 2' Fluro RNA, and chimeric DNA bases

- **Modifications:**
 - 3':** GalNac, Cholesterol, etc.
 - 5':** Palmitate, etc.

- **Purifications:** Desalting, RP-Cartridge, HPLC, PAGE, Dual HPLC

- **Quality Control:** MS and Analytical HPLC

- **Format:** Dried or Customized Buffer



Research and Preclinical



ISO 9001 Certified



2 weeks Turnaround



High Purity, High Yield